





AVNeoTM**-Sizer for Aortic Valve Repair Using Autologous Pericardium-****Product Catalog****Prosthetic Free Aortic Valve Repair
AVNeo**

- Commissural Distance Sizing
- Templates for customized leaflets
- Disposable, easy to use
- 12 sizers for easy sizing of varying aortic valves



For further information visit
avneo.net

For registration of cases, or registry information, visit
avneo.jp

Product Information

AVNeo™ Sizer System



AVNeo™ Sizer System Common Kit (JD-005-CK1)



AVNeo™ Small Sizer (JD-005-SS1)



AVNeo™ Large Sizer (JD-005-LS1)

AVNeo™ Sizer System Common Kit

Sizers 19, 21, 23, 25, 27, 29 (Total 6 sizes)

Templates for molding leaflets

Pericardium Fixation Plates

Petri Dish

0.6% Glutaraldehyde Solution Tray

Leaflet Size Mark Sheet

Qty

Product Number

1

JD-005-CK1

AVNeo™ Small Sizer

13, 15, 17 (3 sizes)

Qty

Product Number

1

JD-005-SS1

AVNeo™ Large Sizer

31, 33, 35 (3 sizes)

Qty

Product Numbers

1

JD-005-LS1

-No exchanges on products used or opened

Official Product Name : OZAKI AVNeo Sizer System

Manufactured and Distributed by :

Japanese Organization for Medical Device Development, Inc.

Nihonbashi Lifescience Building #601

2-3-11 Nihonbashi Honcho Chuo-ku, Tokyo 103-0023

Phone: +81-3-6262-3322 Email: info@jomdd.com

